

INSTRUCTION FOR USE

INSTRUCTION MANUAL

Catalog No: C10920
Components: Wood material.
Packing: 100 pcs.

PROKAMA – WOODEN WEDGE

1- Description

They are wedges made of wood with four different thickness options. Each color represents a different thickness.

2- Usage area / Indications

The wedges provide tooth separation and protect the anatomical tooth and tissue lines of the interproximal region, preventing the matrix from slipping. Concave sections adapt to the interdental contour.

Prokama is compatible with tooth anatomy:

- Triangular edges adapt to the interdental contour.
- Rectangular corners allow for easy insertion and removal of wedges.

3- Contraindications

Do not use in patients with hypersensitivity to its ingredients.

4- Warning:

Care should be taken not to swallow it.

5- Storage conditions

It does not require any special storage conditions.

NOTE:

Our products are produced for use in dentistry. It is not suitable for use other than professional dentists. Since the use of the products is under the control of the user other than us, the user is completely responsible for the applications. We recommend using it in accordance with the information and warnings given in the user manual.

** Disposal: Products should be treated according to their biological hazards and disposed of in accordance with the country's legal regulations and hospital policies.



PROMIDA DIŐ MALZEMELERİ SAN. TİC. LTD. ŐTİ.

75. YIL MAH. KÜÇÜK ORGANİZE SANAYİ 3. CADDE 6A ODUNPAZARI/ESKİŐEHİR

info@promida.com.tr www.promida.com.tr