

INSTRUCTION MANUAL

Catalog No: C10915

Components: Stainless steel material.

Packing: 50 pcs.

PROMATRIX OTO – AUTOMATRIX BAND

1- Description

PROMATRIX OTO is a circumferential matrix system used for posterior teeth. It is produced in 2 different tape thicknesses and 2 different lengths:

- 0.038 mm – 6.3 mm height
- 0.038 mm – 5 mm height

2- Usage Area / Indications

PROMATRIX OTO allows multiple restorations to be completed simultaneously in the same quadrant, eliminating the need for conventional matrix holders. It helps form appropriate proximal contacts more easily and provides sufficient space for restorative material and wedge application. It can be used together with rubber covers or other isolation techniques.

3- Contraindications

Do not use in patients with hypersensitivity to its ingredients.

4- Warning:

Care should be taken not to swallow it.

5- Storage conditions

Store at room temperature. Protect from moisture and extreme heat and cold.

NOTE:

Our products are produced for use in dentistry. It is not suitable for use other than by professional dentists. Since the use of the products is under the control of the user other than us, the user is completely responsible for the applications. We recommend using it in accordance with the information and warnings given in the user manual.

** Disposal: Products should be treated according to their biological hazards and disposed of in accordance with the country's legal regulations and hospital policies.